

LABOR ENDERS

PD Dr. Maren Eggers

c/o Labor Prof. Dr. G. Enders MVZ GbR

Rosenbergstraße 85 • 70193 Stuttgart

Expert Opinion

on the efficacy of

Dr. Heuberger's Desinfektion

against

Poliovirus

PD Dr. Maren Eggers

c/o Labor Prof. Dr. G. Enders MVZ GbR
Rosenbergstraße 85
70193 Stuttgart

LABOR ENDERS

Akkreditiert nach DIN EN ISO/IEC 17025



Labor Prof. Dr. G. Enders MVZ GbR • Rosenbergstraße 85 • 70193 Stuttgart


Dr. Heuberger Medical Products GmbH
Thalenhorststr. 8

28307 Bremen

2020-08-07

The efficacy of **Dr. Heuberger's Desinfektion** against human Poliovirus type 1 strain LSc-2ab was tested in a suspension test according to the European standard DIN EN 14476:2013+A1:2015/prA2:2016. The effectiveness of the disinfectant was evaluated under both clean conditions (0.3 g/l BSA) and dirty conditions (3.0 g/l BSA + 3.0 ml/l erythrocytes) as interfering substance. **Dr. Heuberger's Desinfektion** was tested as a 100 ppm, 1000 ppm and 2000 ppm concentration. The exposure time was 5 minutes.

In conclusion, product Dr. Heuberger's Desinfektion is effective against Poliovirus at room temperature under clean conditions (0.3 g/l BSA) and dirty conditions (3.0 g/l BSA + 3.0 ml/l erythrocytes) as interfering substance with an application time of 5 minutes.



PD Dr. rer. nat. Maren Eggers